

TECNIS

Odyssey™ IOL

Also available in TECNIS™ Toric II

A woman with blonde hair, wearing a light-colored cardigan, is looking down at a tablet device she is holding. The background is a warm, blurred indoor setting, possibly a living room with a fireplace.

**Cataracts
can limit
you, pick
a lens
that won't.**

About TECNIS Odyssey™ and TECNIS Odyssey™ Toric II IOL Physicians considering use of the TECNIS Odyssey™ and TECNIS Odyssey™ Toric II IOLs should refer to the Directions for Use labeling for a complete list of indications and safety information. See reverse for continued Indications, References, and Important Safety Information.

What is a cataract?

Cataracts form as we age and are caused by a buildup of the eye's natural protein. The increase in protein causes a cloudiness in the lens of the eye, which is why people with cataracts may experience blurry or dull vision.¹

Do you have presbyopia?

Presbyopia typically occurs around age 40 and is an age-related loss of the eyes' ability to see objects up close.² Cataract treatment is an opportunity to correct presbyopia at the time of surgery, with a lens like **TECNIS Odyssey™ IOL**, which offers a continuous range of vision*, from up close to far away, day and night.³⁻⁵

Do you have astigmatism?

Astigmatism is when the shape of your eye is more oblong than circular, which can create distorted vision. With **TECNIS Odyssey™ IOL Toric II** astigmatism can be corrected during cataract surgery resulting in less need for glasses.[†] Talk to your doctor about your astigmatism correction options.

Enjoy full range of vision⁵

TECNIS Odyssey™ IOL is a premium intraocular lens used to treat cataracts.

It provides a full range of vision—near, intermediate and far—so you can have freedom to perform different activities throughout the day and night.^{4,5}



Learn more at:
clearvisionforyou.com



Standard Monofocal Intraocular Lens



TECNIS Odyssey™ IOL



Images for illustrative purposes only. Actual results may vary. Please see the back side for Indications and Important Safety Information.

TECNIS

Odyssey™ IOL

Also available in TECNIS™ Toric II

INDICATIONS AND IMPORTANT SAFETY INFORMATION FOR TECNIS ODYSSEY™ AND TECNIS ODYSSEY™ TORIC II IOLS WITH TECNIS SIMPLICITY™ DELIVERY SYSTEM, MODELS DRN00V AND DRT

CAUTION: Federal law restricts this device to sale, distribution and use by or on the order of a physician.

INDICATIONS FOR USE: The TECNIS SIMPLICITY™ Delivery System is used to fold and assist in inserting the TECNIS Odyssey™ IOL, which is indicated for primary implantation for the visual correction of aphakia in adult patients, with less than 1 diopter of pre-existing corneal astigmatism, in whom a cataractous lens has been removed. The TECNIS SIMPLICITY™ Delivery System is used to fold and assist in inserting the TECNIS Odyssey™ Toric II IOLs that are indicated for primary implantation for the visual correction of aphakia and for reduction of refractive astigmatism in adult patients with greater than or equal to 1 diopter of preoperative corneal astigmatism, in whom a cataractous lens has been removed. Compared to an aspheric monofocal lens, the TECNIS Odyssey™ IOLs mitigate the effects of presbyopia by providing improved visual acuity at intermediate and near distances to reduce eyeglass wear, while maintaining comparable distance visual acuity. The lens is intended for capsular bag placement only.

CONTRAINDICATIONS: None.

RISKS: Routine cataract surgery risk, irrespective of lens selection, could be minor, temporary, or affect patients' vision permanently. Rare complications are worsening of vision, bleeding, or infection. Risks related to use of this lens include halo or glare in nighttime or poor visibility conditions. These may be perceived as an annoyance or hindrance, which, on rare occasions, may be significant enough for the patient to request removal of the IOL. Patients should discuss all risks and benefits with their eye doctor before surgery.

PRECAUTIONS: Recent contact lens usage may affect the patient's refraction; therefore, in contact lens wearers, surgeons should establish corneal stability without contact lenses prior to determining IOL power. Prior to surgery, the surgeon must inform prospective patients of the possible risks and benefits associated with the use of this device and provide a copy of the patient information brochure to the patient. Patients with a predicted postoperative astigmatism greater than 1.0 D may not be suitable candidates for implantation with the TECNIS Odyssey™ and TECNIS Odyssey™ Toric II IOLs, as they may not obtain the benefits of reduced spectacle wear or improved intermediate and near vision seen in patients with lower postoperative astigmatism. The safety and effectiveness of the TECNIS Odyssey™ IOL have not been substantiated in patients under the age of 22 or those with preexisting ocular conditions and intraoperative complications.

WARNINGS: Physicians considering lens implantation should weigh the potential risk/benefit ratio. Well-informed patients with well-defined visual needs and preferences should be selected for lens implantation. Patients should be informed of the possibility of visual effects (such as halo or glare around lights), which may be expected in nighttime or poor visibility conditions.

* Continuous 20/32 or better.

† Individual results will vary. Some TECNIS Odyssey™ patients may require spectacles post-surgery.

REFERENCES:

1. National Eye Institute. Cataracts. <https://www.nei.nih.gov/learn-about-eye-health/eye-conditions-and-diseases...> Accessed May 24, 2022.
2. Presbyopia. American Optometric Association Website: <https://www.aoa.org/patients-and-public-eye-and-vision-problems/glossary-of-eye-and-visionconditions/presbyopia>.
3. Data on file. Johnson & Johnson Surgical Vision, Inc. DOF2023CT4007.
4. TECNIS Odyssey™ IOL DFU Z311926E Rev. B 07/2023.
5. Data on file. Johnson & Johnson Surgical Vision, Inc. DOF2023CT4023.